# CLEARFIL™ REPAIR

|           |             |            |  |
|-----------|-------------|------------|--|
| <b>EN</b> | ENGLISH     | English    | DENTAL LIGHT-CURED REPAIR KIT                          |
| <b>FR</b> | FRANÇAIS    | French     | KIT DE RÉPARATION DENTAIRE PHOTOPOLYMÉRISABLE          |
| <b>ES</b> | ESPAÑOL     | Spanish    | KIT DE REPARACIÓN FOTOPOLIMERIZABLE DE USO DENTAL      |
| <b>IT</b> | ITALIANO    | Italian    | KIT DI RIPARAZIONE DENTALE FOTOPOLIMERIZZABILE         |
| <b>NL</b> | NEDERLANDS  | Dutch      | DENTALE LICHTUITHARDENDE REPARATIE SET                 |
| <b>DE</b> | DEUTSCH     | German     | DENTALES lichthärtendes REPARATURSET                   |
| <b>SV</b> | SVENSKA     | Swedish    | DENTAL LJUSHÄRDANDE REPARATIONSSATS                    |
| <b>NO</b> | NORSK       | Norwegian  | DENTALT LYSHERDET REPARASJONSSETT                      |
| <b>FI</b> | SUOMI       | Finnish    | TÄYTTEEN VALOKOVETTEINEN KORJAUSPAKKAUS                |
| <b>DA</b> | DANSK       | Danish     | DENTALT LYSHÆRDENDE REPARATIONS-KIT                    |
| <b>PT</b> | PORTUGUÊS   | Portuguese | KIT DE REPARAÇÃO FOTOPOLIMERIZAVEL PARA USO DENTÁRIA   |
| <b>EL</b> | ΕΛΛΗΝΙΚΑ    | Greek      | ΟΔΟΝΤΙΚΟ ΦΩΤΟΠΟΛΥΜΕΡΙΖΟΜΕΝΟ ΚΙΤ ΕΠΙΔΙΟΡΘΩΣΗΣ           |
| <b>TR</b> | TÜRKÇE      | Turkish    | IŞIKLA SERTLEŞEN DENTAL TAMİR SETİ                     |
| <b>PL</b> | POLSKI      | Polish     | ŚWIATŁOUTWARDZALNY DENTYSTYCZNY ZESTAW NAPRAWCZY       |
| <b>RO</b> | ROMÂNĂ      | Romanian   | KIT DE REPARAȚII DENTARE FOTOPOLIMERIZABIL             |
| <b>HR</b> | HRVATSKI    | Croatian   | DENTALNI SVJETLOSNO POLIMERIZIRANI KOMPLET ZA POPRAVKE |
| <b>HU</b> | MAGYAR      | Hungarian  | FOGÁSZATI, FÉNYRE KEMÉNYEDŐ JAVÍTÓKÉSZLET              |
| <b>SL</b> | SLOVENŠČINA | Slovenian  | DENTALNI SVETLOBNO STRJUJOČI KOMPLET ZA POPRAVILA      |
| <b>CS</b> | ČESKY       | Czech      | DENTÁLNÍ SVĚTLEM TUHNOUCÍ SADA PRO OPRAVY              |
| <b>SK</b> | SLOVENSKY   | Slovak     | DENTÁLNA SVETLOM VYTVRDZOVANÁ SÚPRAVA PRE OPRAVY       |
| <b>BG</b> | БЪЛГАРСКИ   | Bulgarian  | СТОМАТОЛОГИЧЕН ФОТОПОЛИМЕРИЗИРАЩ КОМПЛЕКТ ЗА ПОПРАВКА  |
| <b>UK</b> | УКРАЇНСЬКА  | Ukrainian  | СТОМАТОЛОГІЧНИЙ ФОТОПОЛІМЕРНИЙ РЕСТАВРАЦІЙНИЙ НАБІР    |
| <b>ET</b> | EESTI       | Estonian   | DENTAALNE VALGUSKÕVASTUV PARANDUSKOMPLEKT              |
| <b>LV</b> | LATVIEŠU    | Latvian    | ZOBĀRSTNIECĪBAS GAISMĀ CIETĒJOŠS LABOŠANAS KOMPLEKTS   |
| <b>LT</b> | LIETUVIŠKAI | Lithuanian | ODONTOLOGINIS ŠVIESA KIETINAMAS TAISYMO RINKINYS       |

## Kuraray Noritake Dental Inc.

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

## EMERGO EUROPE

Westervoortsedijk 60,  
6827 AT Arnhem, The Netherlands


## Kuraray Europe GmbH (EU Importer)

Philipp-Reis-Str. 4,  
65795 Hattersheim am Main, Germany  
Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835  
URL: <https://www.kuraraynoritake.eu>

# CLEARFIL™ REPAIR

## I. INTRODUCTION

CLEARFIL REPAIR consists of CLEARFIL SE BOND (PRIMER and BOND), CLEARFIL PORCELAIN BOND ACTIVATOR, ALLOY PRIMER, K-ETCHANT GEL, CLEARFIL ST OPAQUER and accessories. CLEARFIL SE BOND is a light-curing bonding system, and consists of a selfetching primer and a bonding agent. The primer offers simultaneous treatment of both dentin and enamel. CLEARFIL PORCELAIN BOND ACTIVATOR is used in combination with CLEARFIL SE BOND to enhance bonding strength to ceramic or cured composite resin. ALLOY PRIMER is a metal adhesive primer. It is used to enhance the adhesion between dental metals and resin-based materials. K-ETCHANT GEL is an etching gel and consists of 40% phosphoric acid aqueous solution and colloidal silica. CLEARFIL ST OPAQUER is a light-cured restorative composite resin used to mask exposed metal surface and tooth discoloration. The general clinical benefit of CLEARFIL REPAIR is to restore tooth function for the following INDICATIONS FOR USE.

## II. INDICATIONS FOR USE

Intraoral repairs of fractured crowns, bridges, inlays and onlays made of ceramic, hybrid ceramics or composite resin

## III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers or acetone

## IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn whitish when contacted by CLEARFIL SE BOND PRIMER or ALLOY PRIMER due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing. K-ETCHANT GEL may cause inflammation or erosion with its irritancy.

## V. INCOMPATIBILITIES

- Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- Do not use hemostatics containing ferric compounds since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingival due to remaining ferric ions.
- When using hemostatics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherend surface. Failure to do so might weaken the bond strength to the tooth structure.

## VI. PRECAUTIONS

### 1. Safety precautions

- This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
- If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
- If the product comes in contact with human body tissues, take the following actions:
  - <If the product gets in the eye>  
Immediately wash the eye with copious amounts of water and consult a physician.
  - <If the product comes in contact with the skin or the oral mucosa>  
Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol, and rinse with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Do not reuse a disposable brush tip to prevent cross-contamination. Discard it after use.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components.
- If the patient begins to feel sick after inhaling the acetone contained in ALLOY PRIMER, encourage the patient to breathe copious amounts of fresh air.
- Dispose of this product as a medical waste to prevent infection.

### 2. Handling and manipulation precautions

[Common precautions]

- The product must not be used for any purposes other than specified in [II.INDICATIONS FOR USE].
- The use of this product is restricted to dental professionals.
- Use a rubber dam to prevent contamination and to control moisture.

### [CLEARFIL SE BOND]

- Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing. PRIMER and BOND will set to a gel if left under an operating light or natural light (sunlight from windows).
- If there is an adherend surface including uncut enamel, apply K-ETCHANT GEL to the uncut enamel, leave in place for 10 seconds, wash with water, and then dry. Then, apply the product to an adherend surface.
- Clean the cavity sufficiently to prevent poor bonding. If the adherend surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.

### [CLEARFIL PORCELAIN BOND ACTIVATOR]

- If the treated adherend surface is contaminated with saliva or blood, rinse with water, dry, clean the area with an etching agent (e.g. K-ETCHANT GEL), and re-treat again.
- If the product accumulates on the nozzle or the container, remove it by using a cotton pellet or a gauze pad saturated with alcohol.

### [ALLOY PRIMER]

- If the metal adherend surface is contaminated with saliva before or after the application of the product, clean the surface with a cotton pellet moistened with ethanol and apply the primer again.

### [K-ETCHANT GEL]

- If the product adheres on clothes, wash off with water.
- The product may be hardened when it is not in use. In that case, the content returns to liquid as shaking the bottle with the cap.
- Care should be taken to prevent saliva from coming in contact with the etched tooth surface. Any contaminated tooth surface with saliva etc. must be re-etched.

### [CLEARFIL ST OPAQUER]

- Make sure to wipe any remaining paste off the tip of the syringe with a cotton pellet before closing the cap; otherwise the cup might come off.
- Avoid sharing the same paste among different patients to prevent cross-contamination.

### [Dental light-curing unit]

- Do not look directly at the light source. Protective glasses are recommended.
- Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

### 3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- ALLOY PRIMER, K-ETCHANT GEL and CLEARFIL ST OPAQUER must be stored at 2-25°C/36-77°F when not in use. CLEARFIL SE BOND and CLEARFIL PORCELAIN BOND ACTIVATOR must be refrigerated (2-8°C/36-46°F) when not in use, and should be brought to room temperature before using.
- Replace the bottle cap as soon as possible after resin has been dispensed from the bottle.
- ALLOY PRIMER is volatile. The cap must be replaced as soon as possible after use. If the container is left uncapped, the volatile acetone will evaporate.
- The product must be kept away from extreme heat, direct sunlight or a flame. Note that ALLOY PRIMER is flammable. Do not use or store it near open flames.
- The product must be stored in a proper place where only dental practitioners can access.

## VII. COMPONENTS

### 1. Components

Please see the outside of the package for contents and quantity.

- CLEARFIL SE BOND
  - PRIMER
  - BOND
- CLEARFIL PORCELAIN BOND ACTIVATOR
- ALLOY PRIMER
- K-ETCHANT GEL
- CLEARFIL ST OPAQUER
- Accessories
  - Light blocking plate\*
  - Mixing dish\*
  - Disposable blush tips (white)
  - Brush tip handle (white)\*
  - Outer case\*
 The PRIMER and BOND can be easily dispensed using this outer case.
  - \*Consumables

### 2. Ingredients

Principal ingredients

- CLEARFIL SE BOND
  - PRIMER
    - 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (10-30%)
    - 2-Hydroxyethyl methacrylate (HEMA) (10-30%)
    - Hydrophilic aliphatic dimethacrylate (5-10%)
    - dl-Camphorquinone (< 1%)
    - N,N-Diethanol-p-toluidine (< 5%)
    - Water (40-60%)
  - BOND
    - 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (< 10%)
    - Bisphenol A diglycidylmethacrylate (Bis-GMA) (25-45%)
    - 2-Hydroxyethyl methacrylate (HEMA) (20-40%)
    - Hydrophobic aliphatic dimethacrylate (10-30%)
    - dl-Camphorquinone (< 1%)
    - Accelerators (< 3%)
    - N,N-Diethanol-p-toluidine (< 1%)
    - Colloidal silica (5-15%)
- CLEARFIL PORCELAIN BOND ACTIVATOR
  - Hydrophobic aromatic dimethacrylate (30-70%)
  - 3-Methacryloxypropyl trimethoxysilane (30-70%)
- ALLOY PRIMER
  - Acetone (> 90%)
  - 10-Methacryloyloxydecyl dihydrogen phosphate (MDP) (< 0.5%)
  - 6-(4-Vinylbenzyl-N-propyl) amino-1,3,5-triazine-2,4-dithione (< 1%)
- K-ETCHANT GEL
  - Phosphoric acid
  - Water
  - Colloidal Silica
  - Dye
- CLEARFIL ST OPAQUER
  - Silica filler (> 60%)
  - Silanated colloidal silica (< 3%)
  - Colloidal silica (< 3%)
  - Triethyleneglycol dimethacrylate (TEGDMA) (< 10%)
  - Bisphenol A diglycidylmethacrylate (Bis-GMA) (5-25%)
  - dl-Camphorquinone (< 1%)
  - Benzoyl peroxide (< 1%)
  - Titanium dioxide (< 1%)

Units in parentheses are mass %.

## VIII. CLINICAL PROCEDURES

### 1. Preparation of adherend surfaces

Roughen the adherend surfaces using a diamond point, and place a bevel at the marginal area.

### 2. Phosphoric acid treatment of adherend surfaces

Apply K-ETCHANT GEL to the adherend surface (including enamel if present). Leave in place for 5 seconds, then rinse and dry.

### 3. Precious metal surface

When precious metal is used, apply ALLOY PRIMER and dry.

### 4. Silane treatment

- Dispense one drop each of CLEARFIL SE BOND PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the dispensing dish, and mix.

#### [NOTE]

PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

- Apply the mixture to the adherend surface with a disposable brush tip.

- (3) Leave it in place for 5 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces. If the adherend surface extends to the tooth, leave mixture in place for 20 seconds.
- (4) Dry the adherend surface sufficiently by blowing mild air.  
[CAUTION]  
Do not wash treated surface.

#### 5. Bonding

- (1) Dispense the necessary amount of CLEARFIL SE BOND BOND into a well of the dispensing dish immediately before application.
- (2) Apply BOND to the entire cavity wall with a disposable brush tip. Use caution not to allow saliva or exudate to contact the treated surfaces.
- (3) Make a uniform bond film using a gentle air flow.
- (4) Light-cure BOND for 10 seconds with a dental curing light (irradiation wave length: 400-515nm, light intensity: >300mW/cm<sup>2</sup>).

#### 6. Applying the opaquer

Use CLEARFIL ST OPAQUER for metal surface to prevent metal shine through. Light-cure CLEARFIL ST OPAQUER for 40 seconds with a dental curing light (irradiation wave length: 400-515nm, light intensity: >300mW/cm<sup>2</sup>).

#### 7. Place composite resin restorative

Place the composite resin (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY ES Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

#### [WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.


#### [NOTE]

If a serious incident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL", "SE BOND", "CLEARFIL ST" and "CLEARFIL MAJESTY" are registered trademarks or trademarks of KURARAY CO., LTD.