# PermaCem 2.0

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English

#### Product description

PermaCem 2.0 is an all-purpose, self-adhesive, dual-curing composite luting cement for the permanent cementation of crowns, bridges, inlays, onlays and root posts. The use of an additional enamel-dentin bonding agent or of a metal-zirconia primer is not required.

#### Intended purpose

- Definitive cementation of crowns, bridges, inlays, onlays and root posts made of metals (precious and base metals), metal ceramics, composites or all-ceramics (oxide and silicate ceramics, e.g. zirconium oxide or lithium disilicate)
- Definitive cementation of crowns and bridges on implant abutments comprised of titanium or zirconia

#### Limitations of use

Do not use the material if the preparation does not offer sufficient retention (e.g. as in the case of veneers)

#### Indication

Cementation of indirect restorations and prostheses for the treatment of large dental defects, tooth loss, malocclusions and malformations of the teeth.

#### Contraindication

- Do not apply PermaCem 2.0 directly to the exposed dental pulp
- Do not use the material where there are known allergies to any of the components or contact allergies.

#### Patient target group

Persons receiving treatment as part of a dental procedure.

#### Intended user

Dentist

#### Notes for use

- With the Smartmix version the material is applied directly without a dispensing device.
- Ambient light may lead to premature curing of the composite luting cement. Therefore, PermaCem 2.0 should be protected from ambient light after extrusion from the Smartmix Syringe.
- Ensure that excess cement in hard-to-reach areas (proximal, marginal gingival areas) is removed promptly.
- The PermaCem 2.0 curing process can be accelerated. In order to do so, cure PermaCem 2.0 using a suitable curing light.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be at least 600 mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the material.

#### Timing

Working time (including mixing time)*	≤ 60 s
Removal of excess material (without light-curing)	30 to 60 seconds after inserting the restoration or root post
Light-curing time for removal of excess (optional)	1-2 seconds
Maximum intraoral curing time (including mixing time)	420 s
Light-curing time for curing (optional)	At least 20 seconds with a light intensity of $\geq 1000 \text{ mW/cm}^2$ or 40 seconds with a light intensity of $\geq 600 \text{ mW/cm}^2$

\* The specified times are valid at a room temperature of 23 °C and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

#### Using the Smartmix Syringe (for image, see »Handling«)

- Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.
- Attach a new mixing tip [B].

**Note:** In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

**Note:** Make sure that the notches on the mixing tip line up with those on the Smartmix Syringe [C].

3. Lock the mixing tip by turning it 90° clockwise.
4. The material is mixed in the tip when it is expressed and can then be applied directly.

**Note:** When using a Smartmix Syringe for the first time, squeeze out a pea-sized amount of material and discard it!

**Note:** After application, leave the used mixing tip on the Smartmix Syringe as a seal!

## Recommended use

### Cementation of indirect restorations such as crowns, bridges, inlays and onlays

1. Prepare the clean and dry restoration as follows:
  - ▶ Etchable silicate ceramics: Etch the restoration with 5 % hydrofluoric acid in accordance with the instructions from the manufacturer of the etching gel and rinse thoroughly using a water spray. Subsequently dry with air that is free of oil and water, and prepare the surface using a silanizing agent (e.g. Vitique Silane or LuxaBond Universal) in accordance with the manufacturer's instructions.
  - ▶ Oxide ceramics, metals and metal alloys: Pretreat in accordance with the manufacturer's instructions. Unless described otherwise, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, and dry with oil-free and water-free air.
  - ▶ Other materials: Pretreat in accordance with the manufacturer's instructions.
2. If required, clean the cavity or preparation thoroughly using a water spray to remove any residual temporary cement, and then dry. In doing so, ensure that the dentin retains a little residual moisture and does not become too dry.
3. Apply PermaCem 2.0 inside the pretreated restoration.
4. Insert the restoration within 60 seconds after mixing has begun and attach to the dry, prepared tooth using slight pressure.
5. Remove cement residues within 30 to 60 seconds after inserting the restoration with the aid of a microbrush, brush, foam pellet or scaler. Remove interdental excess carefully using dental floss. Alternatively, the cement residues can be briefly (1 - 2 seconds) cured with light and removed with a scaler.
6. Instruct the patient to adopt habitual occlusion.
7. Allow the material to chemically cure completely for 420 seconds after mixing has begun.
8. In the case of light-permeable restorations, additionally cure the material with a suitable light-curing unit for at least 20 seconds at a light intensity of 1000 mW/cm<sup>2</sup> or 40 seconds at 600 mW/cm<sup>2</sup>.

### Cementation of root posts

1. Endodontically prepare the root canal, clean it, and dry it using a paper point. Ensure that the dentin retains a little residual moisture and does not become overly dry.
2. Prepare the clean and dry root post as follows:
  - ▶ Fiber-reinforced root post (e.g. LuxaPost): Clean the root post with ethanol, dry with air that is free of oil and water, and apply silane (e.g. Vitique Silane) where appropriate in accordance with the manufacturer's instructions.
  - ▶ Root posts from other materials: Pretreat in accordance with the manufacturer's instructions.
3. Apply PermaCem 2.0 directly into the prepared canal with the aid of the appropriate tip. In doing so, insert the tip as deeply as possible into the canal.

**Note:** During the application of the composite luting cement, ensure that the tip remains inserted in the material and is only moved upwards with the rising level of the material.

4. Insert the root post within 60 seconds after mixing has begun and cement using slight pressure.
5. Remove cement residues within 30 to 60 seconds after inserting the root post with the aid of a microbrush, brush, foam pellet or scaler.
6. Allow the material to chemically cure completely for 420 seconds after mixing has begun.

**Note:** The material can additionally be light-cured with a suitable light-curing unit for at least 20 seconds at a light intensity of 1000 mW/cm<sup>2</sup> or 40 seconds at 600 mW/cm<sup>2</sup>.

## Interactions

- Residues of methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

## Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

## Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Product performance characteristics

Complies with ISO/TS16506:2017 Class 3 including radiopacity  $\geq 190$  % Al and shade stability.

## Storage/disposal

- Store at 2 – 8 °C/36 – 46 °F in a dry place.
- After application, leave the used mixing tip on the Smartmix Syringe as a seal!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

## Composition

Dental glass, EBPADMA, UDMA, TEDMA, HEMA, MDP, SiO<sub>2</sub>, bis-GMA, additives.  
Inorganic filling materials: approx. 50 vol% (0.02-3 µm).

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