# **Correct Plus™**

Impression Material

# **Correct Plus™ Fast Set**

Impression Material

## **Instructions for Use**

**English**

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Hrvatski

Český

Dansk

Nederlands

Eesti

Suomi

Français

Deutsch

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Magyar

Italiano

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Latviešu

Lietuviu

Norsk

Polski

Português

Românește

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Srpski

Slovenský

Slovensko

Español

Svenska

Türkçe

Manufactured by

**PENTRON**

800-551-0283 / 714-516-7557 / [www.pentron.com](http://www.pentron.com)

## Correct Plus/Correct Plus Fast Set

Correct Plus Hydrophilic Impression Material and Correct Plus Fast Set Impression Material are versatile impression materials designed to accommodate impressions of inlays, crowns and bridges, orthodontic appliances, precision attachments and veneers. The hydrophilic nature of these materials provides outstanding detail in the presence of fluids and exceptional tear strength ensures your detailed impression will remain intact upon removal. Correct Plus Impression Material provides multiple packaging options and various set times to accommodate all of your impression needs. For optimal patient acceptance, Correct Plus Impression Material is available in unflavored and Berry.

### Indications:

- Impression of inlays, crowns and bridges, orthodontic appliances, precision attachments and veneers

### Precautions:

- Due to the fast working time Correct Plus Fast Set is not recommended for full arch impressions. For full arch impressions, use Correct Plus Regular Set Impression Material.

### Storage and Use:

Store and use Correct Plus and Correct Plus Fast Set Impression Materials at ambient temperature.

Refer to outer package for expiration date. Do not use after expiration date.

### Descriptions of viscosities:

The brief descriptions below are intended to help you understand which materials are available and how they might perform.

- **Light Body (not available in Fast Set):** a traditional light viscosity wash material characterized by adequate flowability. Recommended with any tray material in the Correct Plus and Correct Plus Fast Set line.
- **Thick-n-Thin™:** A unique Thixotropic light viscosity wash material that initially stacks on the tooth when syringed then flows into the sulcus to capture critical subgingival detail. Recommended with any tray material in the Correct Plus and Correct Plus Fast Set line.
- **Universal™ Body:** Truly an all-purpose impression material. Use it as a tray material, a light body or a monophase material. A perfect substitute for alginate materials.
- **Tray (also available in Berry flavor; not available in Fast Set):** A heavy viscosity material with a mousse-like consistency and high end

hardness. Ideal for Three-in-One impression trays. Recommended with any wash material in the Correct Plus and Correct Plus Fast Set line.

- **Auto-Mix Putty (also available in Berry Flavor; not available in Fast Set):** A heavy viscosity material that provides the hydraulic drive of a hand-mix putty in a convenient auto-mix cartridge. For use in stock, custom, or triple trays. Recommended with any wash material in the Correct Plus and Correct Plus Fast Set line.
- **Putty (also available in Berry flavor; not available in Fast Set):** An easy to mix, heavy viscosity putty with a high end hardness to resist distortion and provide maximum hydraulic drive of the light body material. Recommended with any wash material in the Correct Plus and Correct Plus Fast Set line.
- **Bite Superfast:** This material was designed exclusively for faster bite registration, with a mousse-like viscosity and high end hardness. When set, its rigidity facilitates trimming.

## Instructions for Use

### Auto Mix Cartridges:

1. Select and/or prepare a rigid tray of sufficient size which provides at least 2-3 mm thickness of space for impression material.
2. Apply Correct VPS™ Tray Adhesive to the tray and allow to air dry for 5 minutes. Adhesive should also be used with perforated trays.
3. Load auto mix cartridges into impression guns, both wash and tray material of choice.
4. Dispose of sealing cap at front end of dual cartridge and extrude 1/4 inch of material to ensure that equal amounts of material are dispensed. Discard. Wipe end of cartridge clean. Attach auto-mix tip. Attach an intraoral tip to the wash material mixing tip.  
Note: Mixing Tips are for single patient use only, to prevent cross-contamination between patients.
5. Dry the field and place cotton rolls.
6. Place your heavy viscosity of choice into impression tray (or see below for Putty mixing instructions, and place putty into tray). **See table below for specific minimum working times.**
7. Remove cotton rolls. Inject syringe mix around tooth preparations. Seat tray.
8. Retain impression in position until set. **See table below for specific maximum intraoral setting times.**
9. Remove the impression from the mouth.
10. Rinse thoroughly under water and disinfect. The impression can be disinfected by soaking it in Glutaraldehyde or other proper disinfectant solutions. Follow the manufacturer's directions.

11. Models may be poured after 15 minutes when using Correct Plus Fast Set Impression Material and after 30 minutes when using Correct Plus Hydrophilic Impression Material.

**Putty:**

1. Using the scoops provided, dispense equal quantities of base and catalyst.
2. Mix with a stiff spatula (or knead with fingers) until a uniform, streak-free color is achieved, usually less than 45 seconds. Optimal performance requires a ratio of 1:1; slight variations in the relative amounts of base or catalyst with neither increase nor decrease work or set times.

**NOTE:** Never mismatch batch sets of base & catalyst, as unpredictable work and/or set times may result. Do not use latex gloves to mix putty. Wash hands thoroughly before mixing putty by hand.

Material Description	Minimum Working Time*	Maximum Intra Oral Setting Time* (from initial mix)
Correct Plus Hydrophilic Impression Material – <b>Auto Mix Cartridges</b>	1 minute	4 minutes and 30 seconds
Correct Plus Hydrophilic Impression Material – <b>Putty</b>	1 minute	4 minutes and 30 seconds
Correct Plus Fast Set Impression Material – <b>Auto Mix Cartridges</b>	45 seconds	2 minutes and 30 seconds
Correct Plus Bite Superfast	15 seconds	1 minute and 15 seconds

\*These times may vary based on storage conditions, temperature, humidity, and age of product.

**Important Notes:**

- Hemostatic agents and retraction cords that contain ferric sulfates may inhibit the set of impression material. Before taking impression, rinse the preparation site thoroughly with water to remove any residual contaminants.
- Due to its rigidity, use caution when using Correct Plus and Correct Plus Fast Set Auto-Mix Putty on periodontally involved teeth or on bridges unless severe undercuts below the pontics are blocked out.

- Correct Plus Hydrophilic Impression Materials and Correct Plus Fast Set Impression Materials should not be intermixed with or used in conjunction with conventional (condensation cured) silicones.
- When handling Putty, certain gloves (such as latex) will inhibit the set. It is suggested the operator mix a small amount of putty to confirm proper setting prior to impression procedure to test for compatibility. Keep jar closed when not in use.
- Correct Plus and Correct Plus Fast Set impressions can be silver or copper-plated in a metalizing bath.

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# **Correct Plus™ Fast Set**

Impression Material

## **LIMITED WARRANTY— LIMITATION OF PENTRON'S LIABILITY**

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product. In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

**PENTRON**

Manufactured by Pentron Clinical

1717 West Collins Ave. Orange, CA 92867 USA  
800-551-0283 / 714-516-7557 / [www.pentron.com](http://www.pentron.com)

**[EC REP]**

European Union Representative  
SpofaDental a.s.

Markova 238

CZ-506 01 Jičín

Phone: +420 493 583 111

DMC1651\_Rev0  
09/2014